


MEDICAL RECOMBINANT HUMAN COLLAGEN CREAM DRESSING

Product of biotechnological fermentation and genetic engineering, belonging to fibrous collagen. It has the structure of a full-spacing protein sequence of 1068, high biological activity, and can adhere to human cells for "error-free" repair.

Macromolecules cover, fill, adhere and repair wounds, plugging the "big holes" in the skin barrier, 100 times larger than animal collagen fragments.

3 Day Control Test After Medical Aesthetics

-21.98%	-20.22%	+35%	93%
Trans-Epidermal Water Loss (TEWL)	Redness	Skin Regeneration Rate	Patient Satisfaction




Package: 1 box* 1.76 oz/50g
Shelf life: 2 Years

Effects

Sensitive Skin/ Hormonal Dermatitis/Dermatitis Eczema (Infant Eczema and Diaper Rash)

Dry Flaky Skin/Redness and Itching/Reducing Wrinkles/Softening Skin

COMPARISON DIAGRAM OF COLLAGEN ATOMIC STRUCTURE




Ingredients

Composed of Recombinant Human Type III Collagen, humectant, emulsifier, water for injection and a small amount of preservatives.



Benefits

- Prevents scars and hyperpigmentation
- Repairs erythema, papules, and couperose
- Restores your skin's healthy glow
- Deeply hydrates overnight for a smooth, refined complexion by morning
- Alleviates allergies, dermatitis, and tightness
- Relieves dryness and itching while keeping skin hydrated all day



LIQUID HUMAN SKIN

New Gen Recombinant
Humanized Collagen

1068



BONNEHEURE


Effective Skincare Product

5 BENEFITS

15S

Stop the Bleeding

- ▶ Stops bleeding and pain
- ▶ Supports the skin Scar-Free Repair (Type III Collagen)
- ▶ Moisturizing and water locking
Brightens skin
(lightens discoloration)



Barrier ▸ Fill ▸ Lead ▸ Reshape

• Barrier

Recombinant human collagen and hyaluronic acid form a mesh structure on the mucosal surface to lock in water and moisturize.

• Fill


When the skin breaks down or atrophies, the cellular gaps in the stratum corneum increase and collagen is able to penetrate the stratum corneum.

• Lead

Replenish the defective Chuan collagen in the dermis to repair the damage while promoting the division and proliferation of fibroblasts after they migrate to the damaged area.

• Reshape

The internal and external synergies of the components work perfectly to bring out the function of active collagen and achieve collagen remodeling.



Collagen Skin Care

Increased Skin Elasticity After Application



- ✓ Robust skin barrier
- ✓ Repairing damaged skin
- ✓ Increased elasticity, tenderness and smoothness

28-day New Collagen-empowered Skin, Validated By Authoritative Institutions

28-75 days = a skin growth cycle

Overall recognition of skin

state improvement: 93.3%


+23.92%

Rehydration and
water-locking
Cuticle water content



-21.98%

Barrier function
Trans epidermal
water loss value



-13.80%

Redness reduction effect
a value in redness zone



28-day New Collagen-empowered Skin, Validated By Authoritative Institutions

28-75 days = a skin growth cycle

Overall recognition of skin

state improvement: 93.3%


+23.92%

Rehydration and
water-locking
Cuticle water content



-21.98%

Barrier function
Trans epidermal
water loss value




-13.80%

Redness reduction effect
a value in redness zone



USE EFFECT



SCOPE OF APPLICATION



Laser plastic surgery

- >Repairing photoelectric damage
- >Rapid relieving redness
- >Reducing blackening after surgery
- >Locking moisture



Surgery trauma

- >Hemostasis and relieving pain
- >Rapidly healing wound
- >Reducing scarring
- >Avoiding pigmentation



Skin diseases

- >Restoring skin function
- >Reducing inflammatory reaction
- >Thickening the cuticle
- >Adjusting water oil balance
- >Repair Acne Damage
- >Avoid acne pits and marks



Current maintenance

- >Supplementing collagen
- >Reducing wrinkle and skin laxity
- >Improving skin texture and complexion
- >Locking skin moisture
- >After Sun Relief
- >Soothes skin after shaving/tattooing

PRODUCT INFORMATION

○ Product Name: | Medical Collagen Functional Dressings for Skin Repair

○ Main Components: | The main ingredient is recombinant human collagen, which is supplemented with moisturizing emollients, co-solvents etc.

○ Packaging Materials: | Medical aluminum tubes

○ Net Weight: | 50g

○ Features: | Restores the barrier function of the skin, and is suitable for treating local skin itching, burning, skin tightness, pain, etc. Caused by dermatitis and eczema. It can repair the skin barrier function and is suitable for dealing with skin dryness and desquamation, erythema, papules, couperose and skin atrophy and thinning caused by damaged barrier function. It can also prevent scar formation and pigmentation.

SERIES PRODUCTS

Project 28-day rate of change	Serum Gel	Repair Spray	Repair Cream	Repair Mask
				
Skin stratum corneum moisture content (Water locking function)	+23.92%	+23.92%	+23.92%	+23.92%
Transepidermal water loss value (Barrier function)	-21.98%	-14.99%	-15.21%	-16.82%
Skin red zone a* value (Redness reduction and sensitization)	-13.80%	-9.90%	-13.09%	-20.22%